


MAMMOMAT Fusion

Where premium meets everyday mammography

siemens-healthineers.com/mammomat-fusion



SIEMENS
Healthineers 



MAMMOMAT Fusion

Where premium meets everyday mammography

Early detection and optimal diagnostic accuracy play a crucial role in mammography screening. Accomplishing this, however, requires a highly standardized workflow and high-quality images. At the same time, mammography screening requires the examination of thousands of women. Financial sustainability is therefore also crucial.

MAMMOMAT Fusion is a premium mammography system made to enhance everyday screening and diagnostics. The refined, automated workflow performs complex tasks at the click of a button, enabling a high patient throughput without compromising on patient experience. The new-generation, high-resolution detector is remarkably robust and dose-efficient, allowing you to enjoy peace of mind, while low operating costs, attractive service packages and high patient throughput help you get your money's worth.

MAMMOMAT Fusion enables the efficiency and accuracy you need by providing a cost-efficient, premium-quality digital mammography system.

Do it in a click

With a selection of premium, easy-to-use features and a refined workflow, you can perform screening and diagnostics quickly and effortlessly.

Enjoy peace of mind

The new, particularly robust detector features trendsetting MaxFill technology, providing high-resolution images to increase your diagnostic confidence at the right dose.

Get your money's worth

With MAMMOMAT Fusion you benefit from an excellent financial performer. Its selection of premium features increases exam speed for high patient throughput. Enjoy low operating costs and attractive dedicated service packages.

Discover how premium technology can improve your everyday screening and diagnostics. Because detail matters.

Do it in a click

Achieve high patient throughput and reduce staff training requirements with a proven workflow and our unique *syngo* user interface.

Ergonomic acquisition workstation

Easily adjust the workstation table¹ to the height of your choice and select your preferred monitor configuration. Enjoy the intuitive operation of one central *syngo*-based workstation.

Complete data integration

Click just once to import patient data from the radiology information system (RIS¹). This saves time and space with no need for a RIS workstation in the examination room. With a huge storage capacity for more than 30,000 images and comprehensive networking, the system easily manages high volumes of image data.

Efficient one-click-to-image

Perform the whole screening procedure with just a click. The particularly intuitive and seamless screening workflow helps you add speed, comfort and efficiency to every examination. By clicking again, the scheduler opens the next patient automatically.

Convenient Single-Touch Positioning

The brief single touch of a button is all you need to move the system between your predefined positions – such as from cranio-caudal to mediolateral-oblique, for time savings you will really notice.

Seamless isocentric rotation


The center of the detector always remains at the same height, so you can move to the next view with just one click, reducing preparation time and adding efficiency.

Multiple X-ray release

A flexible range of release options include the control console button, a handswitch or footswitch, so you can choose the solution that best fits your screening needs. The image is only one click away.

Stereotactic biopsy

Clear-cut procedures and intuitive operations are the key to fast and accurate examinations. With this in mind, we have designed a stereotactic biopsy¹ solution that reduces the necessary workflow steps and guides the operator through the procedure.



Ergonomic acquisition workstation



Convenient Single-Touch Positioning



Stereotactic Biopsy



Enjoy peace of mind

Achieve consistent diagnostic results with the new, particularly robust detector that reliably provides high-resolution images year after year.

New-generation CsI detector

The innovative new-generation CsI detector technology improves spatial resolution and dose efficiency. With its large 23 cm x 30 cm (9" x 12") format, the remarkably robust detector provides the coverage you need and is perfectly suited for challenging situations, such as bus-installed mobile screening units. The system is ready to use quickly, as no warm-up time is required.

Unique MaxFill Technology

MaxFill technology uses the incoming X-rays more efficiently. The effective X-ray-sensitive area of each pixel has been significantly increased with an innovative multi-layer design. This new detector architecture results in a smaller, 83 µm pixel size and higher detective quantum efficiency (DQE) compared to previous technologies, for better image quality at the right dose.

Personalized OpDose®

Depending on breast density and thickness, OpDose determines the exposure parameters required for each breast individually and uses Automatic Exposure Control to acquire the best-possible image. OpDose automatically excludes the pectoral muscle from its calculation.

Personalized OpComp® with SoftSpeed



Thanks to OpComp with SoftSpeed, the compression slows down as soon as the plate reaches the breast. Compression is increased only as long as the breast is soft and pliable, stopping automatically when the optimal level of compression is reached, improving image quality and patient comfort.

Flexible Opview

Enjoy sharper, easily adjustable images using newly refined image processing software. You can select the configuration that best fits your clinical needs from a range of predefined image impressions and flavors.

Calming MoodLight

Getting women to relax during mammograms can be as simple as changing the lighting. The unique MoodLight¹ panel allows you to select specific colors or a series of changing colors.



Get your money's worth

Reduce operating costs while increasing patient throughput with a system focused on financial performance.

The cost-effective choice

A range of workflow-optimizing features, such as one-click-to-image, Single-Touch Positioning and RIS integration as well as a robust and reliable detector technology enable you to screen up to 15 women per hour¹. This adds up to significantly higher throughput and profitability for your screening site.

Wide range of accessories

MAMMOMAT Fusion supports you with magnification, spot imaging, 2D biopsy as well as Stereotactic biopsy² options. A selection of screening and diagnostic compression plates adapt the system to your clinical requirements at any time.

Service tailored to you

Easily customize comprehensive service options from Siemens to your precise requirements. From preventive maintenance to technical phone support and parts coverage, our Protect Plans or Performance Plans are proven service packages that include proactive service and rapid on-site technical response.

Convenient remote service³

In order to meet your performance expectations at any time, we focus on delivering fast response when you need us. Siemens Remote Services as well as our remote updates allow you to optimize your system availability – and your return on investment.



Affordable reading and archiving solution

syngo.plaza⁴ accelerates your mammography reading by loading images quickly and allowing you to define the workflow steps you require. It is an affordable PACS and mammography reading solution that scales to your needs. From a single user workstation to a client-server configuration with floating licenses, it allows multiple users parallel access to all mammography images.



¹ Data on file
² Option
³ Only available to systems that are connected to the Siemens Remote Services platform.


⁴ syngo.plaza is not commercially available in all countries. Due to regulatory reasons it's future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

Why Siemens Healthineers?

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.



Service and exchange

Increasing value by partnering throughout the entire equipment lifecycle

Equipment Maintenance & Monitoring

Reliably servicing your MAMMOMAT Fusion allows you to identify deviations from current norms to maximize equipment availability.

Education Management

Personalized education and training improves your staff's expertise as well as your equipment efficiency.

Fleet Management

A transparent overview of your fleet allows you to manage the performance and maintenance of your Siemens Healthineers equipment, 24/7.

Accessory Solutions

Products from our partner companies complement your use of our equipment in your daily workflow.

Performance Management

An intelligible overview of your radiography performance data helps you make prompt and well-informed decisions.

Asset Management & Planning

Access to innovative medical technology and equipment throughout the entire contract life-time allows you to maximize focus on patient care.

Business Modelling & Financing

Customized business and financial models address your budgetary and enterprise needs enabling you to remain more competitive.

Departmental Layout Optimization

3D-Visualization and digital twin analysis create more efficient workflows and a more enjoyable working environment.

Continuously adding value and caring for your equipment, your staff, your fleet, your workflows, your department, and your entire institution.

MAMMOMAT Fusion

Where premium meets everyday mammography

Calming MoodLight*

An integrated LED panel with a broad range of colors to choose from emits a soothing light. For a more relaxing exam.

Convenient Single-Touch Positioning

moves the system with a single press of a button and allows you to work seamlessly.

Tungsten tube with rhodium filter

saves up to 50% dose compared to Mo/Mo, and enables a high tube load capacity.

Individualized OpComp® with SoftSpeed

automatically stops when the optimal compression for the best possible image quality is reached.

Seamless isocentric rotation

retains the detector at the same height when moving to the next view for faster and easier positioning.

Personalized OpDose® with AEC

individually selects the right dose for each woman.

Ergonomic acquisition workstation

provides access to all patient and exam data on one central *syngo*-based workstation.

Flexible OpView

Selection of five so-called flavors for image impressions. Choose the one for your needs.

Efficient one-click-to-image

starts the exam with a single click and the system is ready for the first exposure.

Stereotactic Biopsy¹

The system intuitively guides you through the entire procedure. Efficient and reliable workflow for optimum diagnostic confidence.



MAMMOMAT Fusion is not commercially available in all countries. Its future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this document are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), it may contain recycled components that adhere to the same quality assurance measures used for factory-new components.

International version.
Not for distribution or use in the U.S.

The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com