

Surgical Tissue Management System Electrosurgical Generator ESG-400



Electrosurgical Generator **ESG-400**

ESG-400 Features and Benefits

Full-range of monopolar and bipolar modes	Perform open, lap and endoscopic surgery procedures	
Tissue-adaptive power output	 Apply the optimal amount of energy Controlled, precise and quick cutting 	
Touch screen	Choose settings easily	
Quick memory function	Find your preferred settings easily	
Instrument recognition	Plug & Play at the Olympus "UNIVERSAL" socket	
New saline modes (for bipolar resection)	Improved ignition performance Continuous vaporization	
Saline detection	Be certain that you are using the right fluid	
Compatibility with Olympus equipment: - USG-400 - FNDOALPHA	THUNDERBEAT instrument OR Integration	
- UHI	Automatic smoke evacuation	

ESG-400 Technical Specifications

Power supply	
Voltage range	100120 V~ / 220240 V~
Frequency	50 / 60 Hz
Maximum input power	1500 VA
Power fuse	10 A (only FST-series from Schurter)
Terminal for potential equalization	Yes

Output

High frequency functions	Monopolar / Bipolar
High frequency	430 kHz ±20%
Maximum high frequency power	320 W
All modes	25% duty cycle
RFCoag (with or without RCAP)	100% duty cycle

Classification

Protection class according to IEC 60601-1	CF, Class I
Classification according to MDD	Ilb
93/42/EEC	

Sockets

MONOPOLAR	2 x 3-pin (Ø 4 mm), Internatonal standard
	1 x 1-pin (Ø 8 mm), Bovie standard
	1 x coaxial (Ø inner 5 mm / Ø outer 9 mm), Erbe standard
BIPOLAR	1 x 2-pin (Ø 4 mm, pin spacing 28.8 mm),
	Internatonal standard
	1 x coaxial (Ø inner 4 mm / Ø outer 8 mm), Erbe standard
UNIVERSAL	1 x 7-pin, Olympus standard
Neutral electrode	2-pin standard, single or split

370 × 465 × 156 mm 12.5 kg



Operation environment

Operation on the orinio	and the same of th
Temperature	+ 10+ 40°C
	(+ 50+ 104°F)
Relative humidity	3085%, non-condensing
Atmospheric pressure	70106 kPa

Output characteristics

		Mode	Maximum output voltage [Vp]	Maximum power level [W]
MONOPOLAR CU	CUT	PureCut	740	300
		BlendCut	1400	200
		PulseCut slow	800	150
		PulseCut fast	800	150
	COAG	SoftCoag	220	200
		PowerCoag	2500	120
		ForcedCoag	2000	120
		SprayCoag	4300	120
BIPOLAR	CUT	BipolarCut	700	100
		SalineCut	700	320
	COAG	BiSoftCoag	220	120
		AutoCoag	220	120
		SalineCoag	230	200
		HardCoag	220	120
		RFCoag	220	50
		FineCoag	220	40

Order information

Article	Article number	Additional information
ESG-400 electrosurgical generator*	WB91051W	
Footswitch double pedal	WB50402W	
Footswitch single pedal	WB50403W	
Power cable, EU, 4.5m, with angled plug	WA95621A	Many European countries, Type E/F
Power cable, US, 4.5m, with angled plug	WA95622A	USA, Canada and other countries, Type B
Power cable, UK, 4.5m, with angled plug	WA95623A	United Kingdom and other countries, Type G
ESG-400 SALINE cable	WA00014A	
Energy cart	TC-E400	
Neutral electrode cable (reusable)	MAJ-814	

^{*} Including one double footswitch.



Width x Depth x Height